

Sani-Cloth Chlor - Information for Use (IFU)

Sani-Cloth Chlor are pre-dosed disposable sporicidal disinfection and cleaning wipes for non-invasive medical devices; non-porous hard surfaces.

Uses: Disinfection for general surfaces in the Medical Area.

Directions: Dispense wipe and open, clean surface in S-motion from top to bottom of area. Ensure the surface remains visibly wet for the required contact time, use another wipe if needed. Allow surface to air dry.

Caution: Avoid contact with eyes, skin and clothing. If irritation occurs discontinue use, if irritation persists seek medical attention. Do not use on invasive medical devices. Follow the medical device manufacturers cleaning and disinfecting guidelines. For Professional Use Only.

PPE: We recommend to using gloves.

Address: PDI Ltd, Aber Park, Flint, Flintshire, CH6 5EX

Storage: Keep product below 23°C. Store upright. Keep out of direct sunlight.

Product Usage: This product is single use. Once a wipe is dispensed from the packaging, it can only be used once and then must be disposed. Overuse and re-use will not provide the efficacy required for correct product use. Do not use if the packaging appears damaged. Do not use any wipes that are dry.

Composition: >5400ppm Sodium Hypochlorite <1% Anionic Surfactant (detergent)

Regulatory Classification: CE marked 93/42/EEC

Product Code	Pack Format	Wipes per Pack	Packs per Case	Wipe Size	NHS Code
XP00297	Canister	50 wipes	6	195x200mm	VJT 228



How to Use the Wipe





1. Check with your internal protocols to ensure that the correct PPE is used for the product, area and equipment you are about to clean and/or disinfect.



Refer to label instructions and 'Use by' date.

2. Check the expiry date.



- - 3. Instructions are available for opening the canisters.





Remove any visible soiling with first wipe before using additional wipes.

4. Wipe use instructions - Remove any visible soiling with the first wipe, then use an additional wipe to disinfect.







5. Remove only one wipe at a time. Open out the wipe.





6. Wipe direction should be 'dirty to clean', top to bottom, taking care not to go over the same area twice to prevent any cross contamination.



Recommended surface area coverage 1m square. Allow disinfectant to dry on surface- Contact time.



Dispose of used wipe into clinical waste.

- 7. Contact Time & Drying Time
 - a. One wipe covers an approximate surface area of one (1) metre square. Do not overuse the wipe, if it becomes dry or soiled discard and use another wipe to complete the area.
 - b. Allow the disinfected area to air dry.
 - c. Dispose of used wipes in the clinical waste bin.
 - d. After use ensure the packaging lid is closed. Once empty, dispose of the packaging in the recycling bin, or according to local protocol



PRODUCT EFFICACY		
EN 16615 4-Field Test (Mechanical Action)	BACTERIA	
Enterococcus hirae	60 sec (c)	
Pseudomonas aeruginosa	60 sec (c)	
Staphylococcus aureus	60 sec (c)	
EN 13727	BACTERIA	
Pseudomonas aeruginosa	30 sec (c)	
Staphylococcus aureus	30 sec (c)	
Enterococcus hirae	30 sec (c)	
EN 16615 4-Field Test (Mechanical Action)	YEAST	
Candida Albicans	1min (c)	
EN 13624	YEAST	
Candida Albicans	30 sec (c)	
Candida Auris	30 sec (c)	
EN 14476	ENVELOPED VIRUS	
Vaccinia virus (HBV, HCV, HIV, H5N1, SARS, Corona)	30 sec (c)	
EN 14476	NON ENVELOPED VIRUS	
Norovirus	60 sec (c)	
Adenovirus	60 sec (c)	
Poliovirus	60 sec (c)	
EN 17126	SPORES	
Bacillus subtilis	15min (c)	
Bacillus cereus	15 min (c)	
Clostridium difficile	15 min (c)	
*EN13704	SPORES	
B subtilis	5min (c)	
EN 16615 4-Field Test (Mechanical Action)	SPORES	
Clostridium difficile	60 sec (c)	
3-Stage Protocol (surface and wipe test)	SPORES	
C difficile	10sec (c) / 5min (d)	

* EN13704 was tested to the new 2018 standard using mature spores

Product Efficacy /Product Claims: the time listed can be followed by letters that indicate the test conditions used during the testing of the product (c) for clean conditions or (d) for dirty conditions.



PRODUCT EFFICACY – EPA Test Results (USA)		
US EPA Guidelines (wipe test)-	BACTERIA	
Pseudomonas aeruginosa ATCC 15442	1min (d)	
Escherichia coli O157:H7 (ATCC 35150)	1min (d)	
Staphylococcus aureus ATCC 6538	1min (d)	
Campylobacter jejuni ATCC 29428	1min (d)	
Klebsiella pneumoniae ATCC 4352	1min (d)	
Legionella pneumophila (ATC 33153)	1min (d)	
Listeria monocytogenes ATCC 19117	1min (d)	
Salmonella Choleraesuis ATCC 10708	1min (d)	
Streptococcus pyogenes ATCC 12344	1min (d)	
Acinetobacter baumannii ATCC 19606	1min (d)	
ESBL Resistant Escherichia coli ATCC BAA - 196	1min (d)	
ESBL Resistant Klebsiella pneumoniae ATCC 700603	1min (d)	
Klebsiella pneumoniae Carbapenem Resistant ATCC BAA - 1705	1min (d)	
Community acquired MRSA (NARSA NRS384) Genotype 300	1min (d)	
Community acquired MRSA (NARSA NRS123) Genotype 400	1min (d)	
Methicillin resistant Staph aureus (MRSA) (ATCC 33591)	1min (d)	
Streptococcus pneumoniae - Penicillium Resistant (ATCC 700677)	1min (d)	
Vancomycin Intermediate Resistant Staphylococcus aureus (VISA	(1)	
HIP 5836)	1min (d)	
Vancomycin Resistant Staphylococcus aureus (VRSA) NARSA VRS1	1min (d)	
Vancomycin Resistant Enterococcus faecalis	1min (d)	
US EPA Guidelines (wipe test)-	FUNGI & YEAST	
Aspergillus niger ATCC 16404	4min (d)	
Candida albicans ATCC 10231 (Yeast)	4min (d)	
Trichophyton mentagrophytes ATCC 9533	4min (d)	
US EPA Guidelines (wipe test)-	ENVELOPED VIRUS	
Hepatitis B (Duck Hep.B virus used as surrogate)	1min (d)	
Hepatitis C (BVD virus used as surrogate)	1min (d)	
HIV - 1 Type 1	1min (d)	
Herpes simplex type 2	1min (d)	
Influenza avian H5N1	1min (d)	
Influenza A Hong Kong	1min (d)	
Human Coronavirus	1min (d)	
Respiratory syncytial virus	1min (d)	
US EPA Guidelines (wipe test)-Southern Research Inst.	NON ENVELOPED VIRUS	



Adenovirus	1min (d)
Poliovirus	1min (d)
Norovirus	1min (d)
Canine parvovirus	1min (d)
Hepatitis A	1min (d)
Rotavirus	1min (d)
Rhinovirus	1min(d)
US EPA Guidelines (wipe test)-	ТВ
Mycobacterium bovis	2min (d)
US EPA Guidelines (wipe test)-	SPORES
C difficile	4min (d)

Product Efficacy /Product Claims: the time listed can be followed by letters that indicate the test conditions used during the testing of the product (c) for clean conditions or (d) for dirty conditions.



Packaging Symbols

Here are some of the most common symbols you may find on PDI product packaging and their description:

Title of Symbol	Description of Symbol	Symbol (Colours may vary)
Do not flush	Do not flush	DO NOT ELLISH
Do not macerate	Do not macerate	DO NOT MACERATE
Consult instructions for use	Indicates the need for the user to consult the Instructions for use, available on the website at <u>www.pdihc.com/global</u>	Ĩ
CE mark	CE mark Medical Device Directive Class Ila	MEDICAL DEVICE
Do not re-use	Indicates a medical device (wipe) that is intended for one use, or for use on a single patient during a single procedure	2
Manufacturer	Indicates the medical device manufacturer, as define in EU directives 90/385/EEC, 93/42EEC and 98/79/EC	
Expressed liquid	Stability and testing completed on liquid expressed from the wipe, rather than bulk liquid.	
ΡΑΟ	Period after opening – the period of time the product can be used once the product has been opened.	1M
Wipe count	How many wipes are in the pack.	50
LOT	Allocated LOT or batch number of the product	LOT
Expiry	The expiry date of the product	



High Density Polyethylene (HDPE)	HDPE 2 is a high density-to-strength ratio and is commonly recycled, and has the number "2" as its resin identification code.	HDPE
-------------------------------------	---	------

Label Elements

Hazard Statements	
H412	Harmful to aquatic life with long lasting effects
Precautionary Statements	
P273	Avoid release to the environment.
P501	Dispose of contents/container in accordance with national regulations.

Authorised Representative in the European Community



NEX MEDICAL ANTISEPTICS S.r.l. Via Per Arluno 37/39 20010 CASOREZZO (MI) Italy



FAQ's

What is contact time and what happens if the surface dries before the stated contact time on a Sani-Cloth product label?

The contact time listed on the **Sani-Cloth** product label is the total amount of time that it takes to inactivate the microorganisms listed on the product label. This time is typically referred to in minutes, and should be communicated to staff members that are utilising the disinfectant. In certain geographies and also in settings where temperature, relative humidity, and air changes may vary, it is possible that the surface may not remain visibly wet for the designated contact time. Current industry guidance requires that the treated environmental surface or equipment remains wet for the contact time stated on product label. Additional wipes may be needed in order to comply with the industry guidance, however the overall contact time does not change.

While the current industry guidance requires the treated environmental surfaces to remain wet for the stated contact time, leading researchers in infection prevention offer an alternate view. In a commentary published in Infection Control and Hospital Epidemiology (March 2018, vol. 39, no. 3, pp 229-331), Dr. W.A. Rutala and Dr. D. J. Weber suggest that contact time and treatment time are mutually exclusive. They suggest that treatment time, irrelevant of wet time, should be followed by healthcare workers for wipes and sprays (except bleach products.) PDI will continue to monitor the science closely and provide their customers with the latest information.

Following cleaning and disinfection the surface should be allowed to fully air dry.

Additional Information

Printed copies of this IFU are available upon request. Email: <u>sales@pdi-hc.co.uk</u> Previous revisions of IFU's are available upon request. Email: <u>sales@pdi-hc.co.uk</u> Document language: English (EN)

Issued 27.07.2020 IFU-10-EN Rev 2