Ready access from the IV pole:









3. Remove individual packets from the bottom to the top of the strip

REMEMBER:

When only one packet is left on strip, remove from IV pole before opening.

***Prevantics**®

3.15% Chlorhexidine Gluconate (CHG) and 70% Isopropyl Alcohol (IPA)



SWAB

• Disinfects needleless access sites prior to use

- Designed for convenience
- Confirmed efficacy at a 5 second scrub time and 5 second dry time



ENVIRONMENT OF CARE PATIENT CARE

INTERVENTIONAL CARE

GENERAL GUIDELINES FOR USE:



70% Isopropyl Alcohol (IPA)

SWAB

1. Look for the words "tear below" (found on back side of packet) and tear downward diagonally to retrieve swab.



2. Prior to access, apply swab to the needleless access site of the device by holding swab between thumb and index finger and scrub using repeated back and forth strokes for 5 seconds.




3. Allow to air dry for 5 seconds. Do not blot or wipe dry.



4. Discard after single use.

5 second scrub time and 5 second dry time.

©2016 PDI Two Nice-Pak Park, Orangeburg, New York 10962 T: 800-444-6725 W: pdihc.com





PDI11168905