

Technical Data Bulletin

Now effective against SARS-CoV-2, the virus that causes COVID-19.

Sani-24[®]



Powered By
MICROBAN



EPA Reg. No. 42182-13-9480 (wipe)

EPA Reg. No. 42182-9-9480 (spray)

Product Description

Sani-24[®] Germicidal Spray and Wipes give you the power with around the clock protection. It is the first, and only, EPA-registered disinfectant with the ability to control HAI-causing microorganisms with Continuously Active Disinfection for up to 24 hours¹.

Chemical Composition

Active Ingredients:

Alkyl dimethyl benzyl ammonium chloride (50% C ₁₄ , 40% C ₁₂ , 10% C ₁₆).....	0.276%
Didecyl dimethyl ammonium chloride.....	0.104%
Octyl decyl dimethyl ammonium chloride.....	0.207%
Diocetyl dimethyl ammonium chloride.....	0.104%
Ethanol.....	68.610%
Other ingredients	30.699%
TOTAL.....	100.000%

Efficacy

Standard Disinfection

Enveloped Viruses – 10 seconds:

<p>Test Method Used:</p> <p>Organic Soil Load:</p> <p>Exposure Time:</p> <p>Results:</p>	<p>Hepatitis B Virus (HBV) (Duck Hepatitis B virus as Surrogate) Virucidal Efficacy of a Disinfectant for Use on Inanimate Environmental Surfaces Whole duck serum (100% duck serum) with an additional 5% fetal bovine serum 10 seconds The results demonstrated complete inactivation of Duck Hepatitis B virus following a 10 second exposure time at 20±1°C (21.0°C), as required by the U.S. EPA</p>
<p>Test Method Used:</p> <p>Organic Soil Load:</p> <p>Exposure time:</p> <p>Results:</p>	<p>Hepatitis C Virus (HCV) (Bovine Viral Diarrhea Virus as Surrogate) [Oregon C24v-genotype 1] Virucidal Efficacy of a Disinfectant for Use on Inanimate Environmental Surfaces 5% horse serum 10 seconds The results demonstrated complete inactivation of Bovine Viral Diarrhea Virus following a 10 second exposure time at 20±1°C (21.0°C), as required by the U.S. EPA</p>
<p>Test Method Used:</p> <p>Organic Soil Load:</p> <p>Exposure time:</p> <p>Results:</p>	<p>Herpes simplex virus type 1 [ATCC VR-733] [F(1)] Virucidal Efficacy of a Disinfectant for Use on Inanimate Environmental Surfaces 5% fetal bovine serum 10 seconds The results indicate complete inactivation of Herpes simplex virus type 1 under these test conditions as required by the U.S. EPA</p>
<p>Test Method Used:</p> <p>Organic Soil Load:</p> <p>Exposure time:</p> <p>Results:</p>	<p>Herpes simplex virus type 2 [ATCC VR-734] [Strain G] Virucidal Efficacy of a Disinfectant for Use on Inanimate Environmental Surfaces 5% fetal bovine serum 10 seconds The results indicate complete inactivation of Herpes simplex virus type 2 under these test conditions as required by the U.S. EPA</p>
<p>Test Method Used:</p> <p>Organic Soil Load:</p> <p>Exposure time:</p> <p>Results:</p>	<p>Severe Acute Respiratory Syndrome-Related Coronavirus 2 (SARS-CoV-2) Virucidal Efficacy of a Disinfectant for Use on Inanimate Environmental Surfaces 5% fetal bovine serum 10 seconds The results indicate complete inactivation of Severe Acute Respiratory Syndrome-Related Coronavirus 2 (SARS-CoV-2) under these test conditions as required by the U.S. EPA</p>
<p>Test Method Used:</p> <p>Organic Soil Load:</p> <p>Exposure time:</p> <p>Results:</p>	<p>Human Coronavirus [ATCC VR-740] [Strain 229E] Virucidal Efficacy of a Disinfectant for Use on Inanimate Environmental Surfaces 5% fetal bovine serum 10 seconds The results indicate complete inactivation of Human Coronavirus under these test conditions as required by the U.S. EPA</p>
<p>Test Method Used:</p> <p>Organic Soil Load:</p> <p>Exposure time:</p> <p>Results:</p>	<p>Human Immunodeficiency virus type 1 (HIV) [Strain HTLV-III_b] Virucidal Efficacy of a Disinfectant for Use on Inanimate Environmental Surfaces 5% fetal bovine serum 10 seconds The results indicate complete inactivation of Human Immunodeficiency virus type 1 (HIV) under these test conditions as required by the U.S. EPA</p>

Avian Influenza A (H3N2) Reassortant virus [ATCC VR-2072]
 [A/Washington/897/80 x A/Mallard/New York/6750/78]
 Test Method Used: Virucidal Efficacy of a Disinfectant for Use on Inanimate Environmental Surfaces
 Organic Soil Load: 5% fetal bovine serum
 Exposure time: 10 seconds
 Results: The results indicate complete inactivation of Avian Influenza A (H3N2) Reassortant virus under these test conditions as required by the U.S. EPA

2009-H1N1 Influenza A virus [(Novel H1N1)] [CDC 2009712192] [A/Mexico/4108/2009]
 Test Method Used: Virucidal Efficacy of a Disinfectant for Use on Inanimate Environmental Surfaces
 Organic Soil Load: 5% fetal bovine serum
 Exposure time: 10 seconds
 Results: The results indicate complete inactivation of 2009-H1N1 Influenza A virus [(Novel H1N1)] under these test conditions as required by the U.S. EPA

Respiratory Syncytial virus (RSV) [ATCC VR-26]
 Test Method Used: Virucidal Efficacy of a Disinfectant for Use on Inanimate Environmental Surfaces
 Organic Soil Load: 5% fetal bovine serum
 Exposure time: 10 seconds
 Results: The results indicate complete inactivation of Respiratory syncytial virus (RSV) under these test conditions as required by the U.S. EPA

Bacteria - 1 minute:

- Acinetobacter baumannii* MDR (Multi-drug resistant) ID#: [ATCC BAA-1605]
- Enterobacter aerogenes* [ATCC 13048]
- Enterobacter aerogenes* MDR (Multi-drug Resistant) ID #: [ATCC 29751]
- Enterococcus faecium* MDR (Multidrug Resistant) ID #: [ATCC 51559]
- Escherichia coli* ESBL (Extended spectrum beta-lactamase) [ATCC BAA-196]
- Escherichia coli* O157:H7 [ATCC 35150]
- Enterococcus faecalis* VRE (Vancomycin resistant enterococcus) [ATCC 51575]
- Klebsiella pneumoniae* CRE (Carbapenem resistant Enterobacteriaceae) [ATCC BAA-2146]
- New Delhi metallo-beta-lactamase-1 (NDM-1) producing Klebsiella pneumoniae (CRE - Carbapenem resistant Enterobacteriaceae) ID #: [ATCC BAA-2146]*
- Pseudomonas aeruginosa* [ATCC 15442]
- Pseudomonas aeruginosa* MBL (Metallo beta-lactamase positive) ID#: [CDC AR-0246/PSA-18]
- Salmonella enterica* [ATCC 10708]
- Staphylococcus aureus* [ATCC 6538]
- Staphylococcus aureus* (Methicillin Resistant) (MRSA) [ATCC 33592]
- Staphylococcus epidermidis* (Methicillin Resistant) (MRSE) [ATCC 51625]
- Staphylococcus aureus* (VISA) (Vancomycin-Intermediate) [HIP5836]
- Staphylococcus aureus* (VRSA) (Vancomycin-Resistant) [HIP11714]

Test Method Used: GLP AOAC Germicidal Spray Products Test, AOAC Germicidal Spray Products as Disinfectants Test modified for towelettes
 Organic Soil Load: 5% fetal bovine serum
 Exposure time: 1 minute
 Incubation: 46-50 hours
 Results: Passed

Fungi-1 minute:

Test Method Used: Fungicidal Germicidal Spray Method, AOAC Fungicidal Germicidal spray products as disinfectants modified for towelettes
 Organic Soil Load: 5% fetal bovine serum
 Exposure time: 1 minute
 Results: No growth observed

Large Non-enveloped

Virus – 2 minutes:

Test Method Used: Virucidal Efficacy of a Disinfectant for Use on Inanimate Environmental Surfaces
 Organic Soil Load: 5% fetal bovine serum
 Exposure time: 2 minutes
 Results: The results indicate complete inactivation of Rotavirus under these test conditions as required by the U.S. EPA.

Fungi – 3 minutes (spray only):

Test Method Used: Fungicidal Germicidal Spray Method
 Organic Soil Load: 5% fetal bovine serum
 Exposure time: 3 minutes
 Results: No growth observed

Fungi – 5 minutes:

Test Method Used: Fungicidal Germicidal Spray Method, AOAC Fungicidal Germicidal spray products as disinfectants modified for towelettes
 Organic Soil Load: 5% fetal bovine serum
 Exposure time: 5 minutes (wipe); 3 minutes (spray)
 Results: No growth observed

Small Non-enveloped

Viruses – 5 minutes:

Test Method Used: Virucidal Efficacy of a Disinfectant for Use on Inanimate Environmental Surfaces
 Organic Soil Load: 5% fetal bovine serum
 Exposure time: 5 minutes
 Results: No growth observed

TB – 5 minutes:

Test Method Used: AOAC Tuberculocidal Activity of Disinfectant Spray Products, AOAC Tuberculocidal Activity of Germicidal Spray Products as Disinfectants Test modified for towelettes
 Organic Soil Load: 5% fetal bovine serum
 Exposure time: 5 minutes (wipe); 3 minutes (spray)
 Results: No growth observed

Continuously Active Disinfection

Bacteria* – 5 minutes:

Acinetobacter baumannii MDR (Multi-drug resistant) ID#: [ATCC BAA-1605]
Enterobacter aerogenes MDR (Multi-drug Resistant) ID#: [ATCC 29751]
Enterobacter aerogenes ATCC 13048
Enterococcus faecalis VRE (Vancomycin resistant enterococcus) ATCC 51575
Enterococcus faecium MDR (Multidrug Resistant) ID#: [ATCC 51559]
 New Delhi metallo-beta-lactamase-1 (NDM-1) producing *Klebsiella pneumoniae* (CRE - Carbapenem resistant *Enterobacteriaceae*) ID #: [ATCC BAA-2146]
Pseudomonas aeruginosa ATCC 15442
Staphylococcus aureus ATCC 6538
Staphylococcus aureus (Methicillin Resistant) (MRSA) ATCC 33592
 Test method used: Modified EPA 01-1A for hospital use claims
 Organic Soil Load: 5% fetal bovine serum
 Dry Contact time: 5 minutes
 Incubation: 46–50 hours
 Results: 99.999% reduction

Hard, Nonporous Non-food Contact Surface Sanitization

Bacteria – 10 seconds:

Enterobacter aerogenes ATCC 13048
Staphylococcus aureus ATCC 6538

Test method:

ASTM E 1153

Organic Soil Load:

5% Fetal Bovine Serum

Incubation:

46–50 hours

Results:

99.9% reduction

Soft Surface Spot Sanitization

Bacteria – 10 seconds:

Enterobacter aerogenes ATCC 13048
Staphylococcus aureus ATCC 6538

Test method:

Modified ASTM E 1153

Organic Soil Load:

5% Fetal Bovine Serum

Incubation:

46–50 hours

Results:

99.9% reduction

TOXICITY

Acute Inhalation

Based on the inhalation test results, **Sani-24** Disinfectant has been classified as Toxicity Category IV for acute inhalation.

Acute Oral Toxicity

Based on the results of this study, **Sani-24** Disinfectant has been classified as Toxicity Category IV for acute oral toxicity.

Acute Eye Irritation

Based on the results of this study, **Sani-24** Disinfectant produced eye irritation that indicates the product would be classified as Toxicity Category II for acute eye irritation.

Acute Dermal Toxicity

Based on the results of this study, **Sani-24** Disinfectant has been classified as Toxicity Category IV for dermal toxicity.

Acute Dermal Irritation

Based on the results of primary skin irritation study, **Sani-24** Disinfectant has been classified as Toxicity Category IV for dermal effects.

Dermal Sensitization

Based upon the sensitization test results, **Sani-24** Disinfectant would not be considered a dermal sensitizing agent.

1. CAD addresses *Acinetobacter baumannii* MDR, *Enterobacter aerogenes*, *Enterobacter aerogenes* MDR, *Enterococcus faecalis* VRE (Vancomycin resistant enterococcus), New Delhi Metallo-beta-lactamase-1 (NDM-1) producing *Klebsiella pneumoniae* (CRE – Carbapenem resistant *Enterobacteriaceae*), *Pseudomonas aeruginosa*, *Staphylococcus aureus*, Methicillin Resistant *Staphylococcus aureus* (MRSA).

* ESKAPE pathogens are *Enterococcus faecium*, *Staphylococcus aureus*, *Klebsiella pneumoniae*, *Acinetobacter baumannii*, *Pseudomonas aeruginosa*, and *Enterobacter* spp

