

# TO BLEACH OR NOT TO BLEACH

PDI, a leader in infection prevention and solutions, understands there are different variables and sensitivities in every healthcare setting; **there is no one wipe which can address every situation.** While not appropriate for facility-wide disinfection, bleach is trusted for disinfecting in critical circumstances.

## PATIENT ROOM

- Mixed patient population
- Constant flux and flow of visitors and healthcare professionals
- Common equipment and surfaces found in healthcare

**Right Product:**  
Super Sani-Cloth®  
Germicidal Disposable Wipe

## OPERATING ROOM

- High oxygen levels
- Airflow sensitivities
- Equipment and hard nonporous surfaces sensitive to alcohol

**Right Product:**  
Sani-Cloth® AF3  
Germicidal Disposable Wipe

## ONCOLOGY

- Highly sensitive patient populations
- High volume of visitors with a desire to keep area clean

**Right Product:**  
Sani-Cloth® AF3  
Germicidal Disposable Wipe

## ISOLATION ROOM

- Confirmed case of *Clostridium difficile* or Norovirus
- Patient on enteric precautions
- Demand for time-tested and trusted disinfectant

**Right Product:**  
Sani-Cloth® Bleach  
Germicidal Disposable Wipe

## EMERGENCY ROOM

- Fast-paced environment, quick room turnover
- Wide range of microorganisms

**Right Product:**  
Super Sani-Cloth®  
Germicidal Disposable Wipe



### To Bleach:

- In cases of extreme outbreaks
- If patients are experiencing symptoms associated with *Clostridium difficile* or Norovirus

Standardizing to one product line enables collaboration between environmental services and nursing staff, which leads to successful infection prevention practices, increased opportunity for cost savings and improved patient outcomes.

**SANI-CLOTH® AF3**  
Germicidal Disposable Wipe  
*Fragrance and Alcohol Free*

Contact Time:  
3 minutes

Effective Against:  
44 Microorganisms

Disinfection Formulation:  
Alcohol Free Quat



**SANI-CLOTH® BLEACH**  
Germicidal Disposable Wipe  
*Trusted for High-Risk Areas*

Contact Time:  
4 minutes

Effective Against:  
50 Microorganisms

Disinfection Formulation:  
1:10 Sodium Hypochlorite



**SUPER SANI-CLOTH®**  
Germicidal Disposable Wipe  
*Fast and Effective*

Contact Time:  
2 minutes

Effective Against:  
30 Microorganisms

Disinfection Formulation:  
Alcohol Quat

