

#1
DISINFECTANT
WIPE IN
HEALTHCARE

Super Sani-Cloth®

GERMICIDAL DISPOSABLE WIPE

Ideal for daily use in fast-paced environments.

- Short **2** minute contact time is effective against enveloped / non-enveloped viruses, bacteria, TB, fungi, multi-drug resistant organisms and bloodborne pathogens²
- Effective against SARS-CoV-2 (Novel Coronavirus), the virus that causes COVID-19
- Designed to be compatible with hard nonporous surfaces and equipment made of stainless steel, plastic, Formica® laminate³, glass and more
- Bactericidal, Tuberculocidal, Virucidal



+PDI®
BE THE DIFFERENCE

Intermediate Level Disinfectant

¹GHX Intelligence Data 2020 Q2

²See Technical Data Bulletin for a complete list of microorganisms


³Formica® is a registered trademark of The Diller Corporation

Providing accessible, compliant, effective disinfection.

Color coded for staff identification



Large compliance icons



Deep well lid provides ample space to store and access next wipe

Easy to read efficacy panel

Point of Care Accessories



PDI® 3-in-1 Universal Sani-Bracket®
for wall/pole mounting of a Sani-Cloth® Canister



Sani-Canister Caddy®
for placement on flat surfaces
(Bilingual)



Sani-TAG® Equipment ID System
to identify what surface disinfectant should be used on a particular item



Educational Wall Charts and Signage
(Bilingual)

	REORDER NO.	WIPE SIZE	CASE PACK	CASE WGT	CASE CUBE	PALLET TI/HI
Super Sani-Cloth® Germicidal Disposable Wipe						
Large Canister	Q55172	6" X 6.75"	160/can, 12 cans/case	22.65 lbs	1.540 ft	10/3
Extra Large Canister	Q86984	7.5" X 15"	65/can, 6 cans/case	14.30 lbs	1.012 ft	10/4
Large Individual Packets	H04082	5" X 8"	50/box, 10 boxes/case	7.72 lbs	0.508 ft	19/5
Extra Large Packets	U87295	11.5" X 11.75"	50/box, 3 boxes/case	8.00 lbs	0.408 ft	24/3

Other formulations available.



Visit pdihc.com for more info.

