


# SKIN + SITE

## Prevent with Prevantics® Antiseptics.

Our swabs feature a **pre-saturated, familiar prep pad** design with a quick **15-second scrub time and 30 second dry time<sup>1</sup>** to help reduce preventable infection rates.

## Prevent with Prevantics® Device Swab.

Our device swab's 3.15% Chlorhexidine Gluconate and 70% Isopropyl Alcohol formula, with **mechanical friction** disinfection, and **fast 5-second scrub and 5-second dry time<sup>2</sup>**, work together to **help lower CLABSI<sup>3</sup> rates.**



**+Prevantics®**

**+PDI®**  
BE THE DIFFERENCE

The **Prevantics** Difference.  
Targeting two critical interventions for preventing CLABSI.



Prevantics Skin Swab



Prevantics Device Swab



	REORDER NO.	CASE PACK	CASE WGT	CASE CUBE	PALLET TI/HI	TREATMENT AREA COVERED
<b>Prevantics® Skin Antiseptic</b>						
Swab 1.0 mL Preinjection indication	B10800	100/box, 10 boxes/case	4.30 lbs	0.27 ft	21/7	2.5 in x 2.5 in (6 cm x 6 cm)
<b>Prevantics® Skin Antiseptic for Kits</b>						
Swab 1.0 mL Preinjection indication	B11400	3000/box, 1 boxes/case	18.93 lbs	1.02 ft	10/4	2.5 in x 2.5 in (6 cm x 6 cm)
<b>Prevantics® Device Swab and Strip</b>						
Swab 1.0 mL	B19600	100/box, 10 boxes/case	4.30 lbs	0.28 ft	21/7	
Swab Strip 1.0 mL	B123ST	160/box, 10 boxes/case	4.68 lbs	0.52 ft	12/7	
Large Swab 1.2 mL	C71200	150/box, 10 boxes/case	7.12 lbs	0.48 ft	12/8	



Learn more



<sup>1</sup>PDI Study Microbac 735-132 July 3, 2014.  
<sup>2</sup>PDI Study BSLI 140304-250 July 23, 2014.  
<sup>3</sup>Central line-associated bloodstream infections.  
© 2024 PDI UPDATE 0124 PDI04238619