

# SKIN + SITE

## Prevent with Prevantics® Antiseptics.

Our Swabsticks and Maxi Swabsticks offer **continued antimicrobial activity** for up to 7 days<sup>1</sup> and feature a low-profile applicator to help staff with **more precise application** to better **reduce preventable infection rates**.

## Prevent with Prevantics® Device Swab.

Our device swab's 3.15% Chlorhexidine Gluconate and 70% Isopropyl Alcohol formula, with **mechanical friction** disinfection, and **fast 5-second scrub and 5-second dry time**<sup>2</sup>, work together to **help lower CLABSI**<sup>3</sup> rates.




**+Prevantics®**

**+PDI®**  
BE THE DIFFERENCE

# The **Prevantics** Difference.

## Targeting two critical interventions for preventing CLABSI<sup>1</sup>s.



REORDER NO.	CASE PACK	CASE WGT	CASE CUBE	PALLET TI/HI	TREATMENT AREA COVERED
-------------	-----------	----------	-----------	--------------	------------------------

### Prevantics® Skin Antiseptics

<b>Swab 1.0 mL</b> Preinjection indication	B10800	100/box, 10 boxes/case	4.30 lbs	0.27 ft	21/7	2.5 in x 2.5 in (6 cm x 6 cm)
<b>Swabstick 1.6 mL</b> Preinjection and preoperative indications	S40750	50/box, 10 boxes/case	7.30 lbs	0.98 ft	15/3	Dry Site: 4 in x 4 in (10 cm x 10 cm) Moist Site: 3 in x 5 in (7.5 cm x 12.5 cm)
<b>Maxi Swabstick 5.1 mL</b> Preinjection and preoperative indications	S41950	30/box, 10 boxes/case	12.90 lbs	1.53 ft	5/5	Dry Site: 8.4 in x 8.4 in (21.3 cm x 21.3 cm) Moist Site: 8.4 in x 8.4 in (21.3 cm x 21.3 cm)

### Prevantics® Skin Antiseptics for Kits

<b>Swab 1.0 mL</b> Preinjection indication	B11400	3000/box, 1 boxes/case	18.93 lbs	1.02 ft	10/4	2.5 in x 2.5 in (6 cm x 6 cm)
<b>Swabstick 1.6 mL</b> Preinjection and preoperative indications	S32450	500/box, 1 boxes/case	5.80 lbs	1.02 ft	10/4	Dry Site: 4 in x 4 in (10 cm x 10 cm) Moist Site: 3 in x 5 in (7.5 cm x 12.5 cm)
<b>Compact Swabstick 1.6 mL</b> Preinjection and preoperative indications	S42850	500/box, 1 boxes/case	6.60 lbs	0.65 ft	10/8	Dry Site: 4 in x 4 in (10 cm x 10 cm) Moist Site: 3 in x 5 in (7.5 cm x 12.5 cm)
<b>Maxi Swabstick 5.1 mL</b> Preinjection and preoperative indications	S23750	300/box, 1 boxes/case	13.20 lbs	1.02 ft	10/4	Dry Site: 8.4 in x 8.4 in (21.3 cm x 21.3 cm) Moist Site: 8.4 in x 8.4 in (21.3 cm x 21.3 cm)

### Prevantics® Device Swab and Strip

<b>Swab 1.0 mL</b>	B19600	100/box, 10 boxes/case	4.30 lbs	0.28 ft	21/7
<b>Swab Strip 1.0 mL</b>	B123ST	160/box, 10 boxes/case	4.68 lbs	0.52 ft	12/7
<b>Large Swab 1.2 mL</b>	C71200	150/box, 10 boxes/case	7.12 lbs	0.48 ft	12/8



**Learn more**



<sup>1</sup>2014 **Prevantics** Clinical Compendium: <https://pdihc.com/resource/prevantics-skin-antiseptics-clinical-compendium-copy/>

<sup>2</sup>PDI Study BSLI 140304-250 July 23, 2014.

<sup>3</sup>Central line-associated bloodstream infections.

© 2024 PDI UPDATE 0124 PD104238620