

#1
BRAND OF
DISINFECTANT
WIPES IN
HEALTHCARE



Sani-Cloth® AF3
GERMICIDAL DISPOSABLE WIPE

Alcohol Free. Fragrance Free. Worry Free.
Protecting your patients has never been easier.

- Short **3** minute contact time is effective against enveloped / non-enveloped viruses, bacteria, TB, fungi, multi-drug resistant organisms and bloodborne pathogens²
- Ideal for use around patients and staff, especially those with respiratory sensitivities
- Alcohol, fragrance, bleach, acid, phenol, acetone, and ammonia free
- Bactericidal, Tuberculocidal, Virucidal



Intermediate Level Disinfectant

¹Market leadership claims based on Q1 2023 Clarivate Data

²See Technical Data Bulletin for a complete list of microorganisms

Available in multiple sizes and formats to meet the various needs of your facility



Refillable pail system:
environmentally conscious and



Easy to read efficacy panel
and bilingual label



Pop-up tissue like
dispensing allows for
one-handed dispensing






160 extra large wipes
for ES convenience

Large compliance icons

80 extra large wipes
ideal for mobile

Point of Care Accessories



**PDI® 3-in-1 Universal
Sani-Bracket®**
for wall/pole mounting of a
Sani-Cloth® Canister



Portable Pack Bracket
for ease of use and accessibility



Sani-Canister Caddy®
color coded for staff identification
and placement on flat surfaces
(Bilingual)



**Educational Wall Charts
and Signage**
(Bilingual)

	REORDER NO.	WIPE SIZE	CASE PACK	CASE WGT	CASE CUBE	PALLET TI/HI
Sani-Cloth® AF3 Germicidal Disposable Wipe						
Large Canister	P13872	6" X 6.75"	160/can, 12 cans/case	26.14 lbs	1.430 ft	10/3
Extra Large Canister	P72584	7.5" X 15"	75/can, 6 cans/case	16.80 lbs	0.936 ft	10/4
Large Individual Packets	H59200	5" X 8"	50/box, 10 boxes/case	7.99 lbs	0.539 ft	20/4
Extra Large Individual Packets	U27500	11.5" X 11.75"	50/box, 3 boxes/case	7.60 lbs	0.500 ft	24/3
Portable Pack	M8063S80	8.2" X 9.8"	80/unit, 9 units/case	18.10 lbs	0.950 ft	10/4
Pail	P1450P	7.5" X 15"	160/pail, 2 pails/case	12.40 lbs	0.953 ft	10/5
Refill (for Pail)	P2450P	7.5" X 15"	160/refill, 2 refills/case	11.10 lbs	0.529 ft	18/4

Other formulations available.



Visit pdihc.com for more info.

400 Chestnut Ridge Road, Woodcliff Lake, NJ 07677
T: 800.999.6423 W: pdihc.com
©2023 PDI UPDATE 0623 PDI01179019