

Full Line Product Catalog



PDI offers a broad range of evidence-based, market-leading *Interventional Care, Environment of Care, and Patient Care* solutions, all designed to help reduce preventable infections, control healthcare costs, and, ultimately, help save lives. It is all part of our passion to Be The Difference[®], every day.



The environment is a primary, and often unexpected, source of infections. Maintain a clean and sanitary environment with our *Environment of Care* solutions, which includes **Sani-Cloth®** disinfectants, the number one brand of disinfecting wipes in healthcare.¹

Sani-Cloth® Prime Germicidal Disposable Wipe & Sani-Prime® Germicidal Spray



Powered by a next generation formulation, protect your patients, staff, and facility from HAIs in just 1 minute.



Sani-HyPerCide™ Germicidal Disposable Wipe & Germicidal Spray



A powerful, ready to use hydrogen peroxide formula designed to address HAI causing microorganisms, including *Clostridioides difficile* without compromising compatibility.



Sani-24® Germicidal Spray



A unique Continuously Active Disinfectant that shields against ESKAPE Pathogens for up to 24 hours.²



Super Sani-Cloth® Germicidal Disposable Wipe



Ideal for daily use in fast-paced environments that require a short contact time and broad coverage of microorganisms.



Sani-Cloth® AF3 Germicidal Disposable Wipe



Alcohol free. Fragrance free. Worry free. Ideal for use around patients and staff, especially those with respiratory sensitivities.



Sani-Cloth® Bleach Germicidal Disposable Wipe



Ideal for disinfecting areas contaminated with *Clostridioides difficile* spores and Norovirus.



These products are on EPA's List N³ for the Emerging Viral Pathogen claim, recommended by the CDC for surface disinfection to help prevent the spread of COVID-19.

¹Market leadership claims based on Q1 2021 Clarivate Data.

²CAD addresses *Acinetobacter baumannii* MDR (Multi-drug resistant), *Enterobacter aerogenes*, *Enterobacter aerogenes* MDR (Multi-drug Resistant), *Enterococcus faecalis* VRE (Vancomycin resistant enterococcus), *Enterococcus faecium* MDR (Multidrug Resistant), New Delhi metallo-beta-lactamase-1 (NDM-1) producing *Klebsiella pneumoniae* (CRE - Carbapenem resistant Enterobacteriaceae), *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Staphylococcus aureus* (Methicillin Resistant) (MRSA)

³<https://www.epa.gov/pesticide-registration/list-n-disinfectants-coronavirus-covid-19>

Sani-Cloth® Plus Germicidal Disposable Cloth



A low-level disinfecting wipe for general surface disinfection.



Easy Screen® Cleaning Wipe



A compatible touchscreen cleaner that works at the speed of technology.⁴

- + Features the power of 70% isopropyl alcohol (IPA)
- + Fast-drying, no residue, anti-fogging, and anti-streaking
- + Rapidly cleans dirt, grime, fingerprints, and smudges⁵
- + Meets AORN's "Guideline for Surgical Attire" to clean cell phones, tablets, and other personal communication or hand-held electronic equipment⁶

UVC Disinfection

Tru-D® Robot



The **Tru-D** robot is a portable UV disinfection system that delivers an automated, measured dose of UVC to consistently disinfect a room, in a single cycle, resulting in the ability to document disinfection results after each and every **Tru-D** room treatment.

- + 20+ validations by independent studies⁷
- + 94% percent decrease in room contamination of epidemiologically important pathogens
- + The **Tru-D** robot operates from one position
- + The **Tru-D** robot's measured dosing capabilities ensure thorough room disinfection
- + Real-time usage-tracking feature provides disinfection documentation to facilities
- + Comprehensive implementation and training for all hospital staff is provided by the **Tru-D** support team

For more information go to tru-d.com or call 800.774.5799

Tru-D® SmartUVC
A PDI SOLUTION

⁴ Consult your manufacturer's instructions for use prior to use.

⁵ Cleaning the screen of your mobile device with a wipe containing alcohol may reduce the life span of any anti-fingerprint coating. Please check your device manufacturer's cleaning instructions.

⁶ AORN Guideline For Surgical Attire, Section 9. Personal Items, Recommendation 9.2, page 12

⁷ <https://tru-d.com/science-papers/>

Different infection risks require different solutions. Our *Interventional Care* solutions help you protect even the most vulnerable patients.

Prevantics® Skin Antiseptics

3.15% (w/v) Chlorhexidine Gluconate (CHG) and 70% (v/v) Isopropyl alcohol (IPA) topical skin antiseptic.

- + Swab is approved for preinjection skin preparation
- + Swabstick and Maxi Swabstick are approved for preoperative and preinjection skin preparation
- + Swabstick and Maxi Swabstick are clinically proven to have continuous antimicrobial activity for up to 7 days⁸
- + Premoistened, ready-to-use applicators contain no glass



Prevantics® Device Swab

3.15% (w/v) CHG/70% (v/v) IPA solution for the disinfection of needleless access sites prior to use and tops of blood culture collection bottles.

- + The only 3.15% CHG/IPA solution to receive 510k clearance from the FDA for disinfecting needleless access sites
- + Strip format hangs conveniently on an IV pole for point-of-care accessibility
- + Prep pad design allows for friction scrub and cleaning of threads rather than just the septum of the needleless connector



Profend® Nasal Decolonization Kit

Easy-to-use, pre-saturated 10% (w/w) PVP-Iodine swabstick – just snap and swab!

- + Assured treatment compliance as product is applied by clinician 1 hour prior to surgery
- + Kills 99.7% of *Staphylococcus aureus* at 1 hour and 99.9% at 12 hours¹¹
- + Ideal for *S. aureus*- and Methicillin-resistant *S. aureus*-colonized patients in the surgical unit, ICU, and other areas of the hospital¹²
- + Neat, dry handle design minimizes mess
- + 4 swabsticks per patient/2 swabsticks per nostril for 15 seconds each



⁸2014 Prevantics® Clinical Compendium.

⁹Dry sites only.

¹⁰Hayden MK, et al. A Randomized Cross-Over Clinical Trial to Compare 3.15% Chlorhexidine/70% Isopropyl Alcohol (CHG) vs. 70% Isopropyl Alcohol Alone (Alcohol) and 5s vs. 15s Scrub for Routine Disinfection of Needleless Connectors (NCs) on Central Venous Catheters (CVCs) in an Adult Medical Intensive Care Unit (ICU), Oral Abstract Presented at 2014 ID Week Conference, October 11, 2014, Philadelphia, PA.

¹¹PDI *in vivo* Study 0113-CTEVO

¹²Stevart D. Antimicrobial-resistant pathogens associated with healthcare associated infections: summary of data reported to the NHSN at the CDC, 2009-2010. *Infect Control Hosp Epidemiol.* 2013;34(1):1-14.

Help your patients and residents remain clean and comfortable with our *Patient Care* solutions.

Sani-Hands® Instant Hand Sanitizing Wipes

Designed for everyone – including staff, patients, residents, and visitors.

- + Contains 70% ethyl alcohol, aloe & vitamin E emollients, and is clinically proven to kill 99.99% of the most common bacteria¹³
- + FDA food code compliant. Safe to use before eating and drinking when used as directed¹⁴
- + Quick drying, leaving no sticky residue



Hygea® Multi-Purpose Washcloths

- + Suitable to use on all areas of the body, including face and neck
- + Premoistened washcloth for full body personal cleansing
- + Gentle cleansing formula contains soothing aloe and vitamin E
- + Hypoallergenic, alcohol-free, and latex-free



Hygea® Flushable Personal Cleansing Cloths

- + Ideal for perineum use and incontinence care
- + Cleans and refreshes better than dry bathroom tissue alone
- + Safe for well-maintained sewers and septic systems¹⁵



¹³Clinical data is representative of Sani-Hands® (NDC #: 10819-5808-1) 65.9% Ethyl Alcohol formulation.

¹⁴Food Code Recommendations of the United States Public Health Service Food and Drug Administration, 2017.

¹⁵Flush only one wipe at a time. Not recommended for use in motor homes or with basement pump systems.

Products kill bugs, but people prevent infections.™ Our point-of-care accessories help put PDI products in the right place at the right time to help healthcare professionals do their job more efficiently and effectively.

All of PDI’s point-of-care accessories are designed to support compliance with facility protocols.

Sani-Bracket®

Easily attach **Sani-Cloth** and **Sani-Hands** canisters to walls or equipment for convenient access. The **Sani-Bracket** is available in Medium, Large/Extra Large, and Large/Extra Large Double Canister.



Sani-Canister Caddy®

The countertop holder for canisters is perfect for nursing stations and carts. It helps display important product and compliance information, such as contact time and precautionary language.



Sani-TAG® Equipment ID System

- + **Sani-TAG** Equipment ID System adheres to surfaces and equipment to help staff identify the specific surface disinfectant that should be used on the particular item.
- + The “Clean Tags” help staff easily identify when equipment has been disinfected and meets regulatory and accreditation standards for identifying clean vs. dirty equipment.
- + All **Sani-TAG** labels remove easily and leave no sticky residue.



PDI Safety Stand

The stand holds one canister of product and has optional attachments that can be affixed directly to the floorstand to hold personal protective equipment.





Reorder Information

REORDER NO. WIPE SIZE CASE PACK CASE WGT CASE CUBE PALLET T/HI

ENVIRONMENT OF CARE



Sani-Cloth Prime Germicidal Disposable Wipe & Sani-Prime Germicidal Spray

Powered by a next generation formulation, protect your patients, staff, and facility from HAIs in just 1 minute.

| | | | | | | |
|---------------------------|--------|----------------|-------------------|-----------|----------|------|
| Large Canister | P25372 | 6" x 6.75" | 12/160s | 16.00 lbs | 1.410 ft | 10/3 |
| X-Large Canister | P24284 | 7.5" x 15" | 6/70s | 27.84 lbs | 1.010 ft | 10/4 |
| Large Individual Packet | H06182 | 5" x 8" | 10/50s | 7.99 lbs | 0.540 ft | 20/4 |
| X-Large Individual Packet | U13195 | 11.5" x 11.75" | 3/50s | 8.65 lbs | 0.500 ft | 24/3 |
| Spray Bottle | X12309 | n/a | 9/32oz per bottle | 19.24 lbs | 0.828 ft | 12/4 |



Sani-24 Germicidal Spray

A unique Continuously Active Disinfectant that shields against ESKAPE pathogens for up to 24 hours.¹⁶

| | | | | | | |
|--------------|--------|-----|-------------------|-----------|----------|------|
| Spray Bottle | X14109 | N/A | 9/32oz per bottle | 20.00 lbs | 1.054 ft | 13/3 |
|--------------|--------|-----|-------------------|-----------|----------|------|



Sani-HyperCide Germicidal Disposable Wipe & Germicidal Spray

A powerful, ready to use hydrogen peroxide formula designed to address HAI causing microorganisms, including *Clostridioides difficile* without compromising compatibility.

| | | | | | | |
|------------------|--------|-----------|-------------------|-----------|----------|------|
| Spray Bottle | X13109 | N/A | 9/32oz per bottle | 22.30 lbs | 1.054 ft | 13/3 |
| Large Canister | P27372 | 6 x 6.75" | 12/160s | 23.6 lbs | 1.360 ft | 10/3 |
| X-Large Canister | P26584 | 7.5 x 15" | 6/65s | 12.5 lbs | 0.862 ft | 10/4 |



Super Sani-Cloth Germicidal Disposable Wipe

Fast, 2-minute contact time allows for a quick room turnover. Ideal for daily use in fast-paced environments that require short contact times and broad coverage of microorganisms.

| | | | | | | |
|---------------------------|--------|----------------|---------|-----------|----------|------|
| Large Canister | Q55172 | 6" x 6.75" | 12/160s | 26.00 lbs | 1.55 ft | 10/3 |
| X-Large Canister | Q86984 | 7.5" x 15" | 6/65s | 15.30 lbs | 1.01 ft | 10/4 |
| Large Individual Packet | H04082 | 5" x 8" | 10/50s | 7.72 lbs | 0.539 ft | 20/4 |
| X-Large Individual Packet | U87295 | 11.5" x 11.75" | 3/50s | 8.00 lbs | 0.50 ft | 24/3 |
| Softpack | A22480 | 8.2" x 9.8" | 9/80s | 15.40 lbs | 0.85 ft | 10/4 |



Sani-Cloth AF3 Germicidal Disposable Wipe

Protecting your staff, patients, and residents has never been easier! Alcohol and fragrance-free formula is ideal for use around those with respiratory sensitivities.

| | | | | | | |
|---------------------------|----------|----------------|---------|-----------|----------|------|
| Large Canister | P13872 | 6" x 6.75" | 12/160s | 26.00 lbs | 1.546 ft | 10/3 |
| X-Large Canister | P63884 | 7.5" x 15" | 6/65s | 15.30 lbs | 1.01 ft | 10/4 |
| Large Individual Packet | H59200 | 5" x 8" | 10/50s | 7.72 lbs | 0.54 ft | 20/4 |
| X-Large Individual Packet | U27500 | 11.5" x 11.75" | 3/50s | 8.00 lbs | 0.50 ft | 24/3 |
| Softpack | M8063S80 | 8.2" x 9.8" | 9/80s | 18.05 lbs | 0.95 ft | 10/4 |
| Pail | P1450P | 7.5" x 15" | 2/160s | 12.40 lbs | 0.95 ft | 10/5 |
| Refill (for Pail) | P2450P | 7.5" x 15" | 2/160s | 11.10 lbs | 0.53 ft | 18/4 |

¹⁶ These 6 pathogens are: *Enterococcus faecium*, *Staphylococcus aureus*, *Klebsiella pneumoniae*, *Acinetobacter baumannii*, *Pseudomonas aeruginosa*, and *Enterobacter* species.



Sani-Cloth Plus Germicidal Disposable Cloth

Ideal for use in alternate care settings, including physician and dental offices.

| | | | | | | |
|------------------|--------|------------|---------|-----------|----------|------|
| Large Canister | Q89072 | 6" x 6.75" | 12/160s | 26.00 lbs | 1.546 ft | 10/3 |
| X-Large Canister | Q85084 | 7.5" x 15" | 6/65s | 15.30 lbs | 1.012 ft | 10/4 |



Sani-Cloth Bleach Germicidal Disposable Wipe

Ideal for disinfecting areas contaminated with *Clostridioides difficile* spores and Norovirus.

| | | | | | | |
|------------------------------|--------|----------------|---------|-----------|----------|------|
| Clinical Size Wipes Canister | P84172 | 6" x 5" | 12/160s | 27.20 lbs | 1.489 ft | 10/3 |
| Large Canister | P54072 | 6" x 10.5" | 12/75s | 26.25 lbs | 1.54 ft | 10/3 |
| X-Large Canister | P25784 | 7.5" x 15" | 6/65s | 20.50 lbs | 1.01 ft | 10/4 |
| Large Individual Packet | H58195 | 5" x 7" | 10/40s | 8.10 lbs | 0.51 ft | 19/5 |
| X-Large Individual Packet | U26595 | 11.5" x 11.75" | 3/40s | 8.20 lbs | 0.50 ft | 24/3 |
| Pail | P7007P | 7.5" x 15" | 2/160s | 16.50 lbs | 0.95 ft | 10/5 |
| Refill (for Pail) | P700RF | 7.5" x 15" | 2/160s | 15.20 lbs | 0.53 ft | 18/3 |



Easy Screen Cleaning Wipe

A compatible touchscreen cleaner that works at the speed of technology. Features the power of 70% IPA and meets AORN's "Guideline for Surgical Attire".

| | | | | | | |
|-----------------|--------|---------|--------|-----------|---------|-----|
| Medium Canister | P03672 | 6" X 9" | 12/70s | 10.64 lbs | .990 ft | 8/5 |
|-----------------|--------|---------|--------|-----------|---------|-----|

INTERVENTIONAL CARE



Profend Nasal Decolonization Kit

Uniquely designed swabsticks pre-saturated with 10% (w/w) Povidone-Iodine (PVP-I) — just snap and swab.

| | | | | | | |
|-------------|--------|-----|---|----------|----------|------|
| Patient Kit | X12048 | n/a | 4 swabs/ patient pack, 12 patient packs/shelf unit, 4 shelf units/case | 2.70 lbs | 0.263 ft | 35/5 |
|-------------|--------|-----|---|----------|----------|------|

REORDER NO. WIPE SIZE CASE PACK CASE WGT CASE CUBE PALLET TI/HI



Prevantics Skin Antiseptics

3.15% (w/v) CHG and 70% (v/v) IPA skin antiseptic.

| | | | | | | |
|------------------------------|--------|-----|---------|-----------|---------|------|
| Swab | B10800 | n/a | 10/100s | 4.30 lbs | 0.27 ft | 21/7 |
| Swabstick | S40750 | n/a | 10/50s | 7.30 lbs | 0.98 ft | 15/3 |
| Maxi Swabstick | S41950 | n/a | 10/30s | 12.90 lbs | 1.53 ft | 5/5 |
| Swab (for kits) | B11400 | n/a | 1/3000s | 18.93 lbs | 1.02 ft | 10/4 |
| Swabstick (for kits) | S32450 | n/a | 1/500s | 5.80 lbs | 1.02 ft | 10/4 |
| Compact Swabstick (for kits) | S42850 | n/a | 1/500s | 6.60 lbs | 0.65 ft | 10/8 |
| Maxi Swabstick (for kits) | S27350 | n/a | 1/300s | 13.20 lbs | 1.02 ft | 10/4 |



Prevantics Device Swab

For disinfecting needleless access sites prior to use. Contains 3.15% (w/v) CHG and 70% (v/v) IPA solution.

| | | | | | | |
|------------|--------|-----|---------|----------|---------|------|
| Swab | B19600 | n/a | 10/100s | 4.30 lbs | 0.28 ft | 21/7 |
| Swab Strip | B123ST | n/a | 10/160s | 4.68 lbs | 0.52 ft | 8/5 |



Alcohol Prep Pads - Sterile

Provides rapid-acting skin antiseptics in a sterile formulation. 70% IPA.

| | | | | | | |
|---------------|--------|-------------|---------|-----------|---------|------|
| Medium Packet | B60307 | 1.1" x 2.6" | 20/200s | 12.40 lbs | 0.84 ft | 8/7 |
| Large Packet | C69900 | 1.7" x 3.5" | 10/100s | 7.20 lbs | 0.48 ft | 12/8 |



PVP Iodine Prep Pads and Swabsticks

Saturated with a 10% (w/w) USP Povidone-Iodine prep solution. A broad-spectrum antiseptic for patient pre-operative skin preparation.

| | | | | | | |
|----------------|--------|-------------|---------|----------|---------|------|
| Medium Packet | B40600 | 1.1" x 2.6" | 10/100s | 4.61 lbs | 0.22 ft | 24/8 |
| Swabstick – 1s | S41350 | n/a | 10/50s | 7.50 lbs | 0.87 ft | 13/4 |
| Swabstick – 3s | S41125 | n/a | 10/25s | 8.70 lbs | 0.87 ft | 13/4 |



PVP Iodine Scrub Swabsticks

Saturated with a 7.5% (w/w) Povidone-Iodine cleansing solution. A broad-spectrum antiseptic for scrubbing skin prior to surgery.

| | | | | | | |
|----------------|--------|-----|--------|----------|---------|------|
| Swabstick – 1s | S48050 | n/a | 10/50s | 7.50 lbs | 0.87 ft | 13/4 |
| Swabstick – 3s | S82125 | n/a | 10/25s | 8.30 lbs | 0.87 ft | 12/4 |



PVP Iodine Duo-Swab® Prep & Scrub Swabsticks

Contains one 7.5% (w/w) Povidone-Iodine cleansing scrub swabstick and one 10% (w/w) Povidone-Iodine antiseptic swabstick in a convenient 2-step package.

| | | | | | | |
|----------------|--------|-----|--------|----------|---------|------|
| Swabstick – 1s | S23125 | n/a | 10/25s | 7.13 lbs | 0.87 ft | 13/4 |
|----------------|--------|-----|--------|----------|---------|------|

PATIENT CARE

**Sani-Hands Instant Hand Sanitizing Wipes**

Contains 70% ethyl alcohol, aloe & vitamin E emollients, and is clinically proven to kill 99.99% of the most common bacteria.¹⁷

| | | | | | | |
|-------------------|--------|-------------|---------|-----------|---------|------|
| Medium Canister | P13472 | 6" x 7.5" | 12/135s | 18.45 lbs | 0.99 ft | 8/5 |
| Large Canister | P15984 | 6" x 7.5" | 6/220s | 15.25 lbs | 0.80 ft | 10/5 |
| Individual Packet | D43600 | 5" x 8" | 10/100s | 9.86 lbs | 0.65 ft | 9/7 |
| Bedside Pack | P71520 | 5.5" x 8.4" | 48/20s | 11.20 lbs | 0.46 ft | 15/6 |

**Hygea Multi-Purpose Washcloths**

Strong, disposable washcloths that are ideal for cleansing and refreshing patients.

| | | | | | | |
|-------------------------|--------|--------------|--------|-----------|---------|------|
| Premium - Solo® Softpak | J14143 | 11.5" x 9.5" | 6/60s | 11.00 lbs | 0.67 ft | 9/7 |
| Resealable Softpak | J22750 | 12" x 8" | 12/48s | 13.30 lbs | 0.67 ft | 10/6 |
| Bulk Individual Packet | U12095 | 11.5" x 10" | 1/400s | 16.70 lbs | 1.56 ft | 6/4 |

**Hygea Flushable Personal Cleansing Cloths**

The ideal flushable pre-moistened washcloth for personal general skin cleansing and incontinence care. Safe for use with sewer and well-maintained septic systems.¹⁸

| | | | | | | |
|--------------|---------|-------------|--------|----------|---------|------|
| Solo Softpak | A500F48 | 5.3" x 6.8" | 12/48s | 6.89 lbs | 0.37 ft | 15/7 |
|--------------|---------|-------------|--------|----------|---------|------|

**Hygea BZK Antiseptic Towelette (alcohol free)**

Multi-purpose antiseptic towelette used for cleansing and refreshing. Contains 0.13% Benzalkonium Chloride.

| | | | | | | |
|-------------------|--------|---------|---------|-----------|---------|-----|
| Individual Packet | D35185 | 7" x 5" | 20/100s | 17.30 lbs | 1.26 ft | 9/3 |
|-------------------|--------|---------|---------|-----------|---------|-----|

**Hygea Obstetrical Towelette**

For personal hygiene and whenever cleansing/refreshing is needed. Saturated with Benzalkonium Chloride, Alcohol 5%, Chlorothymol.

| | | | | | | |
|-------------------|--------|------------|---------|----------|---------|-----|
| Individual Packet | D74800 | 7.75" x 5" | 10/100s | 9.67 lbs | 0.65 ft | 9/7 |
|-------------------|--------|------------|---------|----------|---------|-----|

¹⁷Clinical data is representative of Sani-Hands® (NDC #: 10819-5808-1) 65.9% Ethyl Alcohol formulation.

¹⁸Based on the testing guidelines to determine flushability set by the Association of the Nonwoven Fabrics Industry (INDA). Not for use with basement pump systems. Only flush one wipe at a time.

REORDER NO. WIPE SIZE CASE PACK CASE WGT CASE CUBE PALLET T/H



Hygea Sterile Saline Wipes

Convenient, cost-effective alternative to bottled saline and gauze pads. Delivers a consistent dose of sterile saline with every application.

| | | | | | | |
|-------------------|--------|---------|--------|-----------|---------|------|
| Individual Packet | C22370 | 6" x 4" | 24/24s | 13.37 lbs | 0.49 ft | 23/4 |
|-------------------|--------|---------|--------|-----------|---------|------|



Castile Soap Towelette

Gentle cleansing towelette designed for general cleaning and personal hygiene.

| | | | | | | |
|-------------------|--------|-----------|---------|-----------|---------|-----|
| Individual Packet | D41900 | 7.7" x 5" | 10/100s | 12.00 lbs | 0.65 ft | 9/7 |
|-------------------|--------|-----------|---------|-----------|---------|-----|



Nice 'n Clean® Baby Wipes

Developed for babies' soft and delicate skin and offering a variety of products for every stage of baby care. Available in scented, unscented, and sensitive formulas.

| | | | | | | |
|--|-----------|-------------|--------|-----------|---------|------|
| Scented · Solo® Softpak | M225XT | 7.9" x 6.6" | 12/80s | 12.70 lbs | 0.75 ft | 15/4 |
| Scented · Resealable Softpak | Q34540 | 7.9" x 6.6" | 12/40s | 6.39 lbs | 0.38 ft | 15/8 |
| Unscented · Solo Softpak | M233XT | 7.9" x 6.6" | 12/80s | 12.70 lbs | 0.75ft | 15/4 |
| Unscented · Resealable Softpak | Q70040 | 7.9" x 6.6" | 12/40s | 6.39 lbs | 0.38 ft | 15/8 |
| Newborn Sensitive · Resealable Softpak | Q387S40NY | 7.9" x 6.6" | 12/40s | 6.39 lbs | 0.38 ft | 15/8 |



Sterile Lubricating Jelly

Facilitates the entry of medical and therapeutic devices. Intended for medical purposes. Will not affect surgical instruments, rubber, or plastics. Water soluble, bacteriostatic, latex-free, non-staining, and greaseless.

| | | | | | | |
|----------------------------|--------|-----|---------|-----------|---------|------|
| Individual Packet – 2.7 gm | T00137 | n/a | 12/144s | 13.78 lbs | 1.18 ft | 10/4 |
| Individual Packet – 5.0 gm | T00250 | n/a | 12/72s | 11.40 lbs | 0.90 ft | 12/4 |



Adhesive Tape Remover Pads

Helps to remove adhesive tape from skin and reduce patient discomfort.

| | | | | | | |
|-------------------|--------|------------|---------|----------|---------|------|
| Individual Packet | B16400 | 1.2 x 2.6" | 10/100s | 3.60 lbs | 0.25 ft | 20/6 |
|-------------------|--------|------------|---------|----------|---------|------|



Io-Gone® Iodine Remover Pads

Removes iodine stains from skin and clothing. Quickly dissolves residual iodine color.

| | | | | | | |
|-------------------|--------|-------------|---------|----------|---------|------|
| Individual Packet | B47000 | 1.2" x 2.6" | 10/100s | 4.30 lbs | 0.33 ft | 17/8 |
|-------------------|--------|-------------|---------|----------|---------|------|



See Clear® Lens Cleaning Wipes

Individually wrapped wipe for cleaning protective eyewear, glasses, and loupes.

| | | | | | | |
|-------------------|--------|---------|---------|----------|---------|------|
| Individual Packet | D25431 | 6" x 5" | 12/120s | 9.59 lbs | 0.89 ft | 10/5 |
|-------------------|--------|---------|---------|----------|---------|------|

FOOD SAFETY

**Sani Professional® Cleaning Multi-Surface Wipes**

Ideal for use on hard surfaces, such as dining tables and countertops.

| | | | | | | |
|---------------------------|--------|------------|--------|-----------|------------|------|
| Solo® Softpack (75 count) | A972FW | 9" x 7" | 20/75s | 14.42 lbs | 0.986 ft | 10/3 |
| Solo® Softpack (90 count) | A580FW | 7" x 11.5" | 12/90s | 24.00 lbs | 0.87 ft ft | 10/3 |

POINT-OF-CARE ACCESSORIES

**Sani-Bracket®**

Attach **Sani-Cloth** and **Sani-Hands** canisters to walls or equipment for convenient access; The Medium Bracket, labeled "**Easy Screen**" wipes, will include an optional **Sani-Hands** label. Large/ X-Large Bracket, labeled with "NOT for USE on Skin".

| | | | | | | |
|---|--------|-----|----|-----------|---------|-------|
| PDI Large / Extra Large Canister Sani-Bracket (10 pack) | Q51510 | n/a | 10 | 9.90 lbs | 1.54 ft | 7/4 |
| PDI Large / Extra Large Canister Sani-Bracket (1 pack) | Q51500 | n/a | 1 | 0.88 lbs | 0.11 ft | 56/8 |
| PDI Medium Canister Sani-Bracket (10 pack) | P58510 | n/a | 10 | 9.13 lbs | 1.03 ft | 4/10 |
| PDI Medium Canister Sani-Bracket (1 pack) | P58500 | n/a | 1 | 0.77 lbs | 0.08 ft | 60/10 |
| PDI Large / Extra Large Double Canister Sani-Bracket (10 pack) | Q52510 | n/a | 10 | 19.69 lbs | 3.04 ft | 4/10 |
| PDI Large / Extra Large Double Canister Sani-Bracket (1 pack) | Q52500 | n/a | 1 | 1.65 lbs | 0.22 ft | 60/10 |

**Sani-Canister Caddy®**

Countertop placeholder for canisters. Displays important product and compliance information.

| | | | | | | |
|--|---------|-----|----|----------|---------|-------|
| Sani-Canister Caddy Sani-Hands wipes (10 pack) | P013500 | n/a | 10 | 1.50 lbs | 0.25 ft | 28/6 |
| Sani-Canister Caddy Sani-Hands wipes (5 pack) | P013501 | n/a | 5 | 0.80 lbs | 0.11 ft | 30/10 |
| Sani-Canister Caddy Super Sani-Cloth wipes (10 pack) | P012600 | n/a | 10 | 1.50 lbs | 0.25 ft | 28/6 |
| Sani-Canister Caddy Super Sani-Cloth wipes (5 pack) | P012601 | n/a | 5 | 0.80 lbs | 0.11 ft | 30/10 |
| Sani-Canister Caddy Sani-Cloth AF3 wipes (10 pack) | P015200 | n/a | 10 | 1.50 lbs | 0.25 ft | 28/6 |
| Sani-Canister Caddy Sani-Cloth AF3 wipes (5 pack) | P015201 | n/a | 5 | 0.80 lbs | 0.11 ft | 30/10 |
| Sani-Canister Caddy Sani-Cloth Bleach wipes (10 pack) | P013100 | n/a | 10 | 1.50 lbs | 0.25 ft | 28/6 |
| Sani-Canister Caddy Sani-Cloth Bleach wipes (5 pack) | P013101 | n/a | 5 | 0.80 lbs | 0.11 ft | 30/10 |
| Sani-Canister Caddy Easy Screen wipes | P016101 | n/a | 5 | 0.80 lbs | 0.11 ft | 30/10 |

REORDER NO. WIPE SIZE CASE PACK CASE WGT CASE CUBE PALLET T/H



Softpack Brackets

Attach the softpack to the wall or to mobile equipment for ease of use and accessibility.

| | | | | | | |
|-----------------------------------|---------|-----|----|-----------|---------|-----|
| Softpack Wall Bracket | P039100 | n/a | 10 | 15.50 lbs | 2.34 ft | 6/4 |
| Softpack Mobile Equipment Bracket | P013200 | n/a | 1 | 0.10 lbs | .01 ft | n/a |



Sani-Bracket Deluxe Multi-Canister Floor Stand

Place in entryways and other high traffic areas. Holds up to 5 canisters.

| | | | | | | |
|--|---------|-----|---|-----------|---------|-----|
| Sani-Bracket Deluxe Multi-Canister Floor Stand | P013600 | n/a | 1 | 41.00 lbs | 8.01 ft | 2/2 |
|--|---------|-----|---|-----------|---------|-----|



PDI Safety Stand

Holds one canister of product and has optional attachments that can be affixed directly to the floor stand to hold personal protective equipment.

| | | | | | | |
|------------------|---------|-----|---|-----------|---------|------|
| PDI Safety Stand | P013700 | n/a | 1 | 18.00 lbs | 1.81 ft | 3/7 |
| PPE Station | P014100 | n/a | 1 | 5.50 lbs | 0.64 ft | 11/7 |



Sani-TAG® Equipment ID System

Adheres to equipment to help staff easily identify it as clean. Removes easily and does not leave a sticky residue. Meets regulatory and accreditation standards for identifying clean vs. dirty equipment.

| | | | | | | |
|--|--|-----|--------|----------|---------|------|
| Sani-TAG - Equipment ID System (Clean Tags) | P012200 | n/a | 8/100s | 5.50 lbs | 0.29 ft | 12/5 |
| *Sani-TAG - Disinfect with Sani-Cloth Prime wipes and Sani-Prime spray | *Sani-TAG - Disinfect with Sani-Cloth AF3 wipes | | | | | |
| *Sani-TAG - Disinfect with Sani-24 spray | *Sani-TAG - Disinfect with Sani-Cloth Plus wipes | | | | | |
| *Sani-TAG - Disinfect with Sani-HyPerCide wipes and spray | *Sani-TAG - Disinfect with Sani-Cloth Bleach wipes | | | | | |
| *Sani-TAG - Disinfect with Super Sani-Cloth wipes | *Sani-TAG - Please use Easy Screen wipes | | | | | |

*Request from sales representative

For more information, or to place an order, contact your local PDI sales representative.

©2021 PDI 0921 UPDATE PDI07158138
400 Chestnut Ridge Road, Woodcliff Lake, New Jersey 07677

