

# Sani-24<sup>®</sup> Germicidal Spray

## Product Disclosure

### PRODUCT INFORMATION

Package Size	SKU No.	Country
Spray Bottle (32 oz)	X14109	USA

### INGREDIENT DISCLOSURE

Date of Disclosure: 6/30/21

To learn more:

[California Cleaning Right to Know Act](#)

Ingredients	CAS Numbers	Function	Designated List
Ethanol	64-17-5	Antimicrobial active	None
Water	7732-18-5	Diluent	None
Poly(alkylheteromonocycle)	Not Available	Binder	None
Alkyl C12-16 Dimethylbenzyl Ammonium Chloride	68424-85-1	Antimicrobial active	None
Quaternium-24	32426-11-2	Antimicrobial active	None
Didecyldimonium Chloride	7173-51-5	Antimicrobial active	None
Dimethyldioctylammonium Chloride	5538-94-3	Antimicrobial active	None
Polyalkoxylated Heptamethyltrisiloxane	67674-67-3	Surfactant	None

## SAFETY DATA SHEETS



[Sani-24 Germicidal Spray](#)

## MANUFACTURER

### **PROFESSIONAL DISPOSABLE INTERNATIONAL, INC (PDI)**

400 Chestnut Ridge Road, Woodcliff Lake, NJ 07677

## CUSTOMER SERVICE

(800) 999-6423

Information on date of manufacture can be obtained by calling the toll free number