

SANI-CLOTH® AF3 GERMICIDAL DISPOSABLE WIPE

Product Disclosure

PRODUCT INFORMATION

Package Size	SKU No.	Country
Large Canister (160 count)	P13872	USA
Extra Canister (75 count)	P72584	USA
Large Individual Packets (50 count)	H59200	USA
Extra Large Individual Packets (50 count)	U27500	USA
Softpack (80 count)	M8063S80	USA
Pail (160 count)	P1450P	USA
Refill for pail (160 count)	P2450P	USA

INGREDIENT DISCLOSURE

Date of Disclosure: 1/1/2020

To learn more:

[California Cleaning Right to Know Act](#)

Ingredients	CAS Numbers	Function	Chemicals of Concern
Water	7732-18-5	Solvent	NONE
Dipropylene glycol propyl ether	29911-27-1	Solvent	NONE
Octane-1,2-diol	1117-86-8	Emollient	NONE
Alkyl C12-18 Dimethyl benzyl Ammonium chloride	68391-01-5	Antimicrobial active	NONE

Alkyl C12-18 Dimethyl Ethylbenzyl Ammonium chloride	68956-79-6 85409-23-0	Antimicrobial active	NONE
Octoxynol	9036-19-5	Preservative	NONE
Sodium carbonate	497-19-8	Buffering agent	NONE
Sodium hydroxide	1310-73-2	pH adjuster	OEHHA REL

SAFETY DATA SHEETS



[Sani-Cloth AF3 Germicidal Disposable Wipe](#)

MANUFACTURER

PROFESSIONAL DISPOSABLE INTERNATIONAL, INC (PDI)

400 Chestnut Ridge Road, Woodcliff Lake, NJ 07677

CUSTOMER SERVICE

(800) 999-6423

Information on date of manufacture can be obtained by calling the toll free number