


NEW: Larger treatment area for **Prevantics**[®] Maxi Swabstick!

As of September 1, 2023 the S41950 and S27350 **Prevantics** 5.1mL Maxi Swabsticks will have a larger maximum treatment area, changing from 7" x 7" for dry sites and 3" x 7.5" for moist sites to **8.4" x 8.4" for both dry and moist sites.**

Maximum Coverage Area

Coverage areas below are not to scale



Maxi Swabsticks with updated labeling to reflect the new maximum treatment area will begin to appear over the next several months.

Other important label changes: We've updated our labeling around electrocautery procedures. This change will also appear in the same packaging update as the treatment area size increase.

- + **Previous electrocautery labeling:** Do not use with electrocautery procedures.
- + **NEW electrocautery labeling:** If using an ignition source, allow the solution to completely dry (minimum of 3 minutes on hairless skin; up to 1 hour in hair). Do not blot or wipe away.

Order **Prevantics** Maxi Swabstick today!

- + No glass to break
- + Pre-saturated
- + Unique 2-sided applicator for interdigital areas and skin folds
- + Continued antimicrobial activity for up to 7 days¹



	REORDER NO.	CASE PACK	CASE WGT	CASE CUBE	PALLET TI/HI	TREATMENT AREA COVERED
Prevantics® Skin Antiseptics						
Maxi Swabstick 5.1 mL Preinjection and preoperative indications	S41950	30/box, 10 boxes/case	12.90 lbs	1.53 ft	5/5	Dry Site: 8.4 in x 8.4 in (21.3 cm x 21.3 cm) Moist Site: 8.4 in x 8.4 in (21.3 cm x 21.3 cm)
Prevantics® Skin Antiseptics for Kits						
Maxi Swabstick 5.1 mL Preinjection and preoperative indications	S27350	300/box, 1 boxes/case	13.20 lbs	1.02 ft	10/4	Dry Site: 8.4 in x 8.4 in (21.3 cm x 21.3 cm) Moist Site: 8.4 in x 8.4 in (21.3 cm x 21.3 cm)

⁰On dry sites only.

¹Study: SLM-SC-10, 2004.

²PDI Study Microbac 735-132 July 3, 2014.