

Prevent with Prevantics® Device Large Swab.

Prevantics Device Swabs help improve 'scrub the hub' compliance to better protect against CLABSIs.¹



NEW!

Our device swab's 3.15% Chlorhexidine Gluconate and 70% Isopropyl Alcohol formula, with **mechanical friction** disinfection, and **fast 5-second scrub and 5-second dry time²**, work together to **help lower CLABSI rates.**



+ Prevantics®

+ PDI®
BE THE DIFFERENCE™

NOW AVAILABLE: **NEW Prevantics** Device Large Swab as part of the **Prevantics Device Swab** portfolio.

Features & Benefits

- + 52% larger than current **Prevantics** Device Swab.
- + Conveniently designed for larger vascular access sites.
- + Reduced risk of accidental contact with needleless access site.
- + Quick 5 second scrub time and 5 second dry time helps with staff compliance.
- + Achieves a $>4.0 \log_{10}$ reduction against 7 tested clinically relevant microorganisms known to cause HAIs:^{2,3}




- *Candida albicans*
- *Escherichia coli*
- Methicillin Resistant *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- *Candida parapsilosis*
- *Staphylococcus aureus*
- *Staphylococcus epidermidis*

How to Use


Familiar, easy to use prep pad format.

Step 1




Tear down diagonally.

Step 2



5 second scrub time.

Step 3



5 second dry time.

	REORDER NO.	CASE PACK	CASE WGT	CASE CUBE	PALLET TI/HI
Prevantics® Device Large Swab					
Large Swab 1.2 mL	C71200	150/box, 10 boxes/case	7.12 lbs	0.48 ft	12/8

Learn more



¹Central line-associated bloodstream infections.

²PDI Study BSLI 140304-250 July 23, 2014.

³CDC HAI Prevalence Survey. Magill, S.S., Edwards, J.R., Bamberg, W., et al. Multistate Point Prevalence Survey of Healthcare-Associated Infections. N Engl J Med 2014;370:1198-208.

©2024 PDI 400 Chestnut Ridge Road Woodcliff Lake, NJ 07677

T: 800-999-6423 W: pdihc.com PDI01248957

