

Environment of Care Product Portfolio



Sani-HP1™
GERMICIDAL WIPE

DISINFECTS IN **1** MINUTE

HYDROGEN PEROXIDE FORMULATION

KILLS **23** MICROORGANISMS

Pending Full Regulatory Approval¹

ONE MINUTE, ONE STEP, ONE SMART CHOICE FOR DAILY DISINFECTION—INNOVATIVE HYDROGEN PEROXIDE FORMULATION WITH EXCELLENT SURFACE COMPATIBILITY AND NO PPE REQUIRED²

Sani-Cloth® Prime
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN **1** MINUTE

DISINFECTANT QUAT/ALCOHOL FORMULATION

KILLS **57** MICROORGANISMS

1 MINUTE OVERALL CONTACT TIME, INTERMEDIATE LEVEL DISINFECTANT FOR ALL YOUR HIGH-TOUCH, HOUSE-WIDE SURFACE DISINFECTION NEEDS

Sani-HyPerCide®
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN **1** MINUTE

EFFECTIVE AGAINST C. DIFF³

KILLS **24** MICROORGANISMS

STANDARDIZE YOUR DISINFECTION PROTOCOLS WITH THIS POWERFUL H₂O₂ FORMULA THAT KILLS *CLOSTRIDIODES DIFFICILE* WITHOUT COMPROMISING COMPATIBILITY

Super Sani-Cloth®
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN **2** MINUTES

#1 DISINFECTANT WIPE IN HEALTHCARE³

DISINFECTANT QUAT/ALCOHOL FORMULATION

KILLS **34** MICROORGANISMS

IDEAL FOR DAILY USE IN FAST-PACED ENVIRONMENTS THAT REQUIRE SHORT CONTACT TIMES AND BROAD COVERAGE OF MICROORGANISMS

Sani-Cloth® AF3
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN **3** MINUTES

DISINFECTANT QUAT ALCOHOL FREE FORMULATION

KILLS **45** MICROORGANISMS

PROTECTING YOUR PATIENTS HAS NEVER BEEN EASIER. FRAGRANCE FREE – IDEAL FOR USE AROUND PATIENTS AND STAFF, ESPECIALLY THOSE WITH RESPIRATORY SENSITIVITIES

Sani-Cloth® Plus
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN **3** MINUTES

DISINFECTANT QUAT FORMULATION

KILLS **18** MICROORGANISMS

SIMPLE, RELIABLE PROTECTION—ONE-STEP DISINFECTION IDEAL FOR PHYSICIAN OFFICES, DENTAL CLINICS, AND OTHER ALTERNATE CARE SETTINGS

Sani-Cloth® Bleach
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN **4** MINUTES

DISINFECTANT 1:10 BLEACH DILUTION FORMULATION

KILLS **52** MICROORGANISMS

IDEAL FOR DISINFECTING AREAS CONTAMINATED WITH *CLOSTRIDIODES DIFFICILE* SPORES, *CANDIDA AURIS*, AND NOROVIRUS

Easy Screen®
CLEANING WIPE

IDEAL FOR ANY EQUIPMENT REQUIRING 70% IPA

CLEANING 70% ISOPROPYL ALCOHOL FORMULATION

A COMPATIBLE 70% IPA TOUCHSCREEN CLEANER THAT WORKS AT THE SPEED OF TECHNOLOGY⁵

*Kills SARS-CoV-2 (Coronavirus), the virus that causes COVID-19, on hard, nonporous surfaces excluding Easy Screen® wipes.

*Effective against *Candida auris* (*C. auris*).

¹State regulatory approvals are in progress. Product availability may vary—contact your PDI rep for availability in your state.; ²For all PDI products, according to the label instruction, the use of gloves or other PPE is not required to handle the product, particularly in non-clinical settings. You should, however, wear PPE as appropriate per your facility protocol.; ³Market leadership claims based on Q4 2024 Clarivate Data. Based on units sold to acute and non-acute healthcare settings.; ⁴Sani-Cloth® Plus Germicidal Disposable Wipe compares to competitor Quat/Low-Alcohol disinfectant products.; ⁵Cleaning the screen of your mobile device with a wipe containing alcohol may reduce the lifespan of any anti-fingerprint coating. Please check your device manufacturer's cleaning instructions before using Easy Screen® wipes and all other PDI disinfectant products. Easy Screen® wipes are to be used on energized electronic equipment only.

Sani-HP1™

GERMICIDAL WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Medium Canister	P62572	6" x 5"	275/can 12 cans/case
Large Canister	P42672	6" x 9"	165/can 12 cans/case
Extra-Large Canister	P83184	7.5" x 15"	95/can 6 cans/case
Bulk	P0125P	11" x 12"	300/bag 2 bags/case

Sani-Cloth® Prime

GERMICIDAL DISPOSABLE WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Large Canister	P25372	6" x 6.75"	160/can 12 cans/case
Extra-Large Canister	P24284	7.5" x 15"	70/can 6 cans/case
Large Individual Packets	H06182	5" x 8"	50/box 10 boxes/case
X-Large Individual Packets	U13195	11.5" x 11.75"	50/box 3 boxes/case

Sani-HyPerCide®

GERMICIDAL DISPOSABLE WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Large Canister	P27372	6" x 6.75"	160/can 12 cans/case
Extra-Large Canister	P26584	7.5" x 15"	65/can 6 cans/case

Super Sani-Cloth®

GERMICIDAL DISPOSABLE WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Large Canister	Q55172	6" x 6.75"	160/can 12 cans/case
Extra-Large Canister	P86984	7.5" x 15"	75/can 6 cans/case
Large Individual Packets	H04082	5" x 8"	50/box boxes/case
X-Large Individual Packets	U87295	11.5" x 11.75"	50/box 3 boxes/case
Softpack	A22480	8.2" x 9.8"	80/pack 9 packs/case

Sani-Cloth® AF3

GERMICIDAL DISPOSABLE WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Large Canister	P13872	6" x 6.75"	160/can 12 cans/case
Extra-Large Canister	P72584	7.5" x 15"	75/can 6 cans/case
Large Individual Packets	H59200	5" x 8"	50/box boxes/case
X-Large Individual Packets	U27500	11.5" x 11.75"	50/box 3 boxes/case
Pail	P1450P	7.5" x 15"	160/pail 2 pails/case
Refill (for Pail)	P2450P	7.5" x 15"	160/refill 2 refills/case

Sani-Cloth® Plus

GERMICIDAL DISPOSABLE WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Large Canister	Q89072	6" x 6.75"	160/can 12 cans/case
Extra-Large Canister	Q85084	7.5" x 15"	65/can 6 cans/case

Sani-Cloth® Bleach

GERMICIDAL DISPOSABLE WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Clinical Size Wipes Canister	P84172	6" x 5"	160/can 12 cans/case
Large Canister	P54072	6" x 10.5"	75/can 12 cans/case
Extra-Large Canister	P25784	7.5" x 15"	65/can 6 cans/case
Large Individual Packets	H58195	5" x 7"	40/box 10 boxes/case
X-Large Individual Packets	U26595	11.5" x 11.75"	40/box 3 boxes/case
Pail	P7007P	7.5" x 15"	160/pail 2 pails/case
Refill (for Pail)	P700RF	7.5" x 15"	160/refill 2 refills/case

Easy Screen®

CLEANING WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Medium Canister	P03672	6" x 9"	70/can 12 cans/case
Packet	H01050	5" x 8"	50/box 10 boxes/case

Educational and point-of-care accessories.
Developed to assist you in facility protocol compliance.



PDI Sani-Bracket®
for wall/pole mounting of a canister



Sani-Canister Caddy®
for placement on flat surfaces (bilingual)



Sani-TAG® Equipment ID System
for identification of equipment to ensure the proper wipe is used on a particular piece of equipment



Educational Wall Charts and signage (bilingual)



Scan here to contact your local PDI representative.