

# SANI-CLOTH® AF3 GERMICIDAL DISPOSABLE WIPE

## Product Disclosure

### PRODUCT INFORMATION

| Package Size                              | SKU No.  | Country |
|-------------------------------------------|----------|---------|
| Large Canister (160 count)                | P13872   | USA     |
| Extra Canister (75 count)                 | P72584   | USA     |
| Large Individual Packets (50 count)       | H59200   | USA     |
| Extra Large Individual Packets (50 count) | U27500   | USA     |
| Softpack (80 count)                       | M8063S80 | USA     |
| Pail (160 count)                          | P1450P   | USA     |
| Refill for pail (160 count)               | P2450P   | USA     |
| Bulk packaging (250 count)                | P938250  | USA     |

### INGREDIENT DISCLOSURE

**Date of Disclosure:** 1/1/2020

To learn more:

[California Cleaning Right to Know Act](#)

| Ingredients                                          | CAS Numbers               | Function             | Designated List(s)                    |
|------------------------------------------------------|---------------------------|----------------------|---------------------------------------|
| Water                                                | 7732-18-5                 | Solvent              | NONE                                  |
| Dipropylene glycol propyl ether                      | 29911-27-1                | Solvent              | NONE                                  |
| Caprylyl Glycol                                      | 1117-86-8                 | Wetting agent        | NONE                                  |
| Alkyl C12-18<br>Dimethyl benzyl<br>ammonium chloride | 53516-76-0/<br>68391-01-5 | Antimicrobial active | <a href="#">CA Priority Chemicals</a> |

|                                                           |                           |                      |                                             |
|-----------------------------------------------------------|---------------------------|----------------------|---------------------------------------------|
| Alkyl C12-18<br>Dimethyl ethylbenzyl<br>ammonium chloride | 68956-79-6/<br>85409-23-0 | Antimicrobial active | <a href="#">CA Priority<br/>Chemicals</a>   |
| Octoxynol                                                 | 9036-19-5                 | Surfactant           | <a href="#">EU Endocrine<br/>Disruptors</a> |
| Polyethylene glycol                                       | 25322-68-3                | Surfactant           | NONE                                        |
| Sodium carbonate                                          | 497-19-8                  | Buffering agent      | NONE                                        |
| Sodium hydroxide                                          | 1310-73-2                 | pH adjuster          | <a href="#">OEHA REL</a>                    |

## SAFETY DATA SHEETS



[Sani-Cloth AF3 Germicidal Disposable Wipe](#)

## MANUFACTURER

### **PROFESSIONAL DISPOSABLE INTERNATIONAL, INC (PDI)**

400 Chestnut Ridge Road, Woodcliff Lake, NJ 07677

## CUSTOMER SERVICE

(800) 999-6423

Information on date of manufacture can be obtained by calling the toll-free number