

Sani-Cloth® Bleach Germicidal Disposable Wipe

Product Disclosure

PRODUCT INFORMATION

Package Size	SKU No.	Country
Medium Wipe Canister (160 count)	P84172	USA
Large Canister (75 count)	P54072	USA
Extra Large Canister (65 count)	P25784	USA
Large Individual Packets (40 count)	H58195	USA
Extra Large Individual Packets (40 count)	U26595	USA
Pail (160 count)	P7007P	USA
Refill for pail (160 count)	P700RF	USA
Bulk packaging (220 count)	P753220	USA

INGREDIENT DISCLOSURE

Date of Disclosure: 1/1/2020

To learn more:

[California Cleaning Right to Know Act](#)

Ingredients	CAS Numbers	Function	Designated List(s)
Water	7732-18-5	Solvent	NONE
Trisodium phosphate dodecahydrate	10101-89-0	Chelating agent	NONE
Sodium hypochlorite	7681-52-9	Antimicrobial active	NONE
Sodium hydroxide	1310-73-2	pH adjuster	OEHHA REL
Phosphonobutanetricarboxylic Acid	37971-36-1	Stabilizer	NONE
Ethoxylated phosphate ester	68130-47-2	Surfactant	NONE

SAFETY DATA SHEETS



Sani-Cloth Bleach Germicidal Disposable Wipe

MANUFACTURER

PROFESSIONAL DISPOSABLE INTERNATIONAL, INC (PDI)

400 Chestnut Ridge Road, Woodcliff Lake, NJ 07677

CUSTOMER SERVICE

(800) 999-6423

Information on date of manufacture can be obtained by calling the toll-free number