

# Environment of Care Product Portfolio



**Sani-HP1™**  
GERMICIDAL WIPE

DISINFECTS IN **1** MINUTE

HYDROGEN PEROXIDE FORMULATION

KILLS **58** MICROORGANISMS



ONE MINUTE, ONE STEP, ONE SMART CHOICE FOR DAILY DISINFECTION—INNOVATIVE HYDROGEN PEROXIDE FORMULATION WITH EXCELLENT SURFACE COMPATIBILITY AND NO PPE REQUIRED<sup>1</sup>

**Sani-Cloth® Prime**  
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN **1** MINUTE

DISINFECTANT QUAT/ALCOHOL FORMULATION

KILLS **57** MICROORGANISMS



1 MINUTE OVERALL CONTACT TIME, INTERMEDIATE LEVEL DISINFECTANT FOR ALL YOUR HIGH-TOUCH, HOUSE-WIDE SURFACE DISINFECTION NEEDS

**Sani-HyPerCide®**  
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN **1** MINUTE

EFFECTIVE AGAINST **5** MINUTES

KILLS **24** MICROORGANISMS



STANDARDIZE YOUR DISINFECTION PROTOCOLS WITH THIS POWERFUL H<sub>2</sub>O<sub>2</sub> FORMULA THAT KILLS *CLOSTRIDIODES DIFFICILE* WITHOUT COMPROMISING COMPATIBILITY

**Super Sani-Cloth®**  
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN **2** MINUTES

#1 DISINFECTANT WIPE IN HEALTHCARE<sup>2</sup>

DISINFECTANT QUAT/ALCOHOL FORMULATION

KILLS **34** MICROORGANISMS



IDEAL FOR DAILY USE IN FAST-PACED ENVIRONMENTS THAT REQUIRE SHORT CONTACT TIMES AND BROAD COVERAGE OF MICROORGANISMS

**Sani-Cloth® AF3**  
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN **3** MINUTES

DISINFECTANT QUAT ALCOHOL FREE FORMULATION

KILLS **45** MICROORGANISMS



PROTECTING YOUR PATIENTS HAS NEVER BEEN EASIER. FRAGRANCE FREE - IDEAL FOR USE AROUND PATIENTS AND STAFF, ESPECIALLY THOSE WITH RESPIRATORY SENSITIVITIES

**Sani-Cloth® Plus**  
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN **3** MINUTES

DISINFECTANT QUAT FORMULATION<sup>3</sup>

KILLS **18** MICROORGANISMS



SIMPLE, RELIABLE PROTECTION—ONE-STEP DISINFECTION IDEAL FOR PHYSICIAN OFFICES, DENTAL CLINICS, AND OTHER ALTERNATE CARE SETTINGS

**Sani-Cloth® Bleach**  
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN **4** MINUTES

DISINFECTANT 1:10 BLEACH DILUTION FORMULATION

KILLS **52** MICROORGANISMS



IDEAL FOR DISINFECTING AREAS CONTAMINATED WITH *CLOSTRIDIODES DIFFICILE* SPORES, *CANDIDA AURIS*, AND NOROVIRUS

**Easy Screen®**  
CLEANING WIPE

IDEAL FOR ANY EQUIPMENT REQUIRING 70% IPA

CLEANING 70% ISOPROPYL ALCOHOL FORMULATION



A COMPATIBLE 70% IPA TOUCHSCREEN CLEANER THAT WORKS AT THE SPEED OF TECHNOLOGY<sup>4</sup>

\*Kills SARS-CoV-2 (Coronavirus), the virus that causes COVID-19, on hard, nonporous surfaces excluding Easy Screen® wipes.

\*Effective against *Candida auris* (*C. auris*).

<sup>1</sup>For all PDI products, according to the label instruction, the use of gloves or other PPE is not required to handle the product, particularly in non-clinical settings. You should, however, wear PPE as appropriate per your facility protocol.; <sup>2</sup>Market leadership claims based on Q3 2025 Clarivate data units and dollars sold to healthcare facilities.; <sup>3</sup>Sani-Cloth® Plus Germicidal Disposable Wipe compares to competitor Quat/Low-Alcohol disinfectant products.; <sup>4</sup>Cleaning the screen of your mobile device with a wipe containing alcohol may reduce the lifespan of any anti-fingerprint coating. Please check your device manufacturer's cleaning instructions before using Easy Screen® wipes and all other PDI disinfectant products. Easy Screen® wipes are to be used on energized electronic equipment only.



## Sani-HP1™

GERMICIDAL WIPE

FORMAT	REORDER NO.	WIPE SIZE	CASE PACK
Medium Wipe	P62572	6" x 5"	12 canisters, 275 each
Large Wipe	P42672	6" x 9"	12 canisters, 165 each
X-Large Wipe	P83184	7.5" x 15"	6 canisters, 95 each
Bulk	P0125P	11" x 12"	2 bulk bags, 300 each

## Sani-Cloth® Prime

GERMICIDAL DISPOSABLE WIPE

FORMAT	REORDER NO.	WIPE SIZE	CASE PACK
Large Wipe	P25372	6" x 6.75"	12 canisters, 160 each
X-Large Wipe	P24284	7.5" x 15"	6 canisters, 70 each
Large Packet	H06182	5" x 8"	10 boxes, 50 each
X-Large Packet	U13195	11.5" x 11.75"	3 boxes, 50 each

## Sani-HyPerCide®

GERMICIDAL DISPOSABLE WIPE

FORMAT	REORDER NO.	WIPE SIZE	CASE PACK
Large Wipe	P27372	6" x 6.75"	12 canisters, 160 each
X-Large Wipe	P26584	7.5" x 15"	6 canisters, 65 each

## Super Sani-Cloth®

GERMICIDAL DISPOSABLE WIPE

FORMAT	REORDER NO.	WIPE SIZE	CASE PACK
Large Wipe	Q55172	6" x 6.75"	12 canisters, 160 each
X-Large Wipe	P86984	7.5" x 15"	6 canisters, 75 each
Large Packet	H04082	5" x 8"	10 boxes, 50 each
X-Large Packet	U87295	11.5" x 11.75"	3 boxes, 50 each
Softpack	A22480	8.2" x 9.8"	9 packs, 80 each

## Sani-Cloth® AF3

GERMICIDAL DISPOSABLE WIPE

FORMAT	REORDER NO.	WIPE SIZE	CASE PACK
Large Wipe	P13872	6" x 6.75"	12 canisters, 160 each
X-Large Wipe	P72584	7.5" x 15"	6 canisters, 75 each
Large Packet	H59200	5" x 8"	10 boxes, 50 each
X-Large Packet	U27500	11.5" x 11.75"	3 boxes, 50 each
Bulk	P938250	11" x 12"	2 bulk bags, 250 each

## Sani-Cloth® Plus

GERMICIDAL DISPOSABLE WIPE

FORMAT	REORDER NO.	WIPE SIZE	CASE PACK
Large Wipe	Q89072	6" x 6.75"	12 canisters, 160 each
X-Large Wipe	Q85084	7.5" x 15"	6 canisters, 65 each

## Sani-Cloth® Bleach

GERMICIDAL DISPOSABLE WIPE

FORMAT	REORDER NO.	WIPE SIZE	CASE PACK
Clinical Wipe	P84172	6" x 5"	12 canisters, 160 each
Large Wipe	P54072	6" x 10.5"	12 canisters, 75 each
X-Large Wipe	P25784	7.5" x 15"	6 canisters, 65 each
Large Packet	H58195	5" x 7"	10 boxes, 40 each
X-Large Packet	U26595	11.5" x 11.75"	3 boxes, 40 each
Pail	P7007P	7.5" x 15"	2 pails, 160 each
Refill (for Pail)	P700RF	7.5" x 15"	2 refills, 160 each

## Easy Screen®

CLEANING WIPE

FORMAT	REORDER NO.	WIPE SIZE	CASE PACK
Medium Wipe	P03672	6" x 9"	12 canisters, 70 each
Packet	H01050	5" x 8"	10 boxes, 50 each

### Educational and point-of-care accessories.

Developed to assist you in facility protocol compliance.



PDI Sani-Bracket® for wall/pole mounting of a canister



Sani-Canister Caddy® for placement on flat surfaces (bilingual)



Sani-TAG® Equipment ID System for identification of equipment to ensure the proper wipe is used on a particular piece of equipment



Educational Wall Charts and signage (bilingual)



Scan here to contact your local PDI representative.